



Contribution ID : 57

Type : **Poster**

## Recent developments on Echidna

*Monday, 4 November 2024 17:35 (20)*

Echidna has just installed an improved monochromator, providing both increased intensity for common wavelengths and convenient access to a new wavelength, 2.05Å at the standard 140° degree takeoff angle. The new monochromator is characterised in this poster, together with insights into Echidna usage and exploration of performance limits, including: how little sample is needed? What is the detection limit for magnetic moments? How absorbing can my sample be?

### Topics

Biological Systems and Soft Matter

**Primary author(s) :** AVDEEV, Max (Australian Nuclear Science and Technology Organisation, Australian Centre for Neutron Scattering); Dr HESTER, James (The Australian Nuclear Science and Technology Organisation ); Dr WANG, Chin-Wei (ANSTO); MAYNARD-CASELY, Helen (Australian Nuclear Science and Technology Organisation); PETERSON, Vanessa (ANSTO)

**Presenter(s) :** Dr HESTER, James (The Australian Nuclear Science and Technology Organisation )

**Session Classification :** Posters